



Methead Replacement: 1st Ray



Methead Replacement 1st MPJ

Implant for the replacement of the metatarsal head in the first MPJ.

General Specifications:

Articular surface is an spherical. Stem is a cannulated screw with deep thread.

Material: Cobalt Chrome. All screws are hand sharpened.

Thread D sizes 1 & 2	6.5 mm
Thread D sizes 3 to 5	7.5 mm
Cannulation all sizes	1.6 mm

Style No.	Base OD	Shaft length
MH14-16	14	11.3
MH16-16	16	12.5
MH18-16	18	13.8
MH20-16	20	14.8
MH22-16	22	16.0



System components:


- K100-16S 1.6x100 mm trocar K-wire, single ended
- X240 Met-Head CHI tray
- H100 Handle for shaft driver
- H200 Handle for shaft driver, w/swivel
- X630-16-13 Shaft met head driver



Methead Replacement: 1st Ray



Methead Replacement 1st MPJ



- | | |
|-------------|----------------------------|
| Z720-16 | K-wire holder |
| M2SZ-xx | Doubl-ended sizer xx:14-22 |
| X110 | Driver pin protector |
| X240 | Met Head, CHI tray |
| X620-16 | CHI Driver shaft |
| X630-16-13 | Met Head Shaft Driver |
| Z400-45-150 | Cannulated drill |



X630-16-13 methead driver