

This implants for the replacement of the hemi or metatarsal head in lesser MPJ. Although the metatarsal head implanted are packaged together with hemi implants for the replace of the base of the proximal phalanx, under no circumstance should both a hemi and a metatarsal head be used in the same joint.

**General Specifications:**

Articular surface is spherical. Stem is a cannulated screw with deep thread.

All screws are hand sharpened.

All mini hemis are available in titanium only and have 1.1 mm cannulation.

Style No.	Base D	Material
MH09-11	9	CoCr
MH10-11	10	CoCr
MH11-11	11	CoCr
MH12-11	12	CoCr
MH13-11	13	CoCr
MCHI-08T	8	Titanium
MCHI-09T	9	Titanium
MCHI-10T	10	Titanium
MCHI-11T	11	Titanium
MCHI-12T	12	Titanium





**System components:**

- K100-11S 1.1x100 mm K-wire, single ended
- H100 Handle for shaft driver
- H200 Handle for shaft driver, w/swivel
- MHSZ-09 9 mm sizer
- MHSZ-10 10 mm sizer
- MHSZ-11 11 mm sizer
- MHSZ-12 12 mm sizer
- MHSZ-13 13 mm sizer
- X300 Mini Met-Head, CHI tray
- X600-11-78 Mini hemi driver shaft
- X610-11-75 Mini met head driver shaft
- Z720-11 K-wire holder